

## Attaching an ear mould to a hearing aid

If the hospital has posted your child a new ear mould you will receive a package like this:



Your child's name should be printed here.

**Please check it is the correct name.**

1. Wash your hands.



2. Remove the old hearing aid mould by gripping the ear hook in one hand and pulling on the tubing with the other:



3. Remove the new ear mould from the packaging and take the tubing and push it onto the ear hook of the hearing aid.



4. Check in your child's ear. If it is too long, trim it using clean scissors, approximately **3mm** at a time, until you have found the correct length.



5. Repeat until the hearing aid is positioned correctly on the ear.

